

## Application Note

**Product Name** : KAPA HiFi HotStart DNA polymerase  
**Manufacturer** : KAPA BIOSYSTEMS  
**Application** : Fidelity verification of KAPA HiFi HotStart DNA polymerase by Human TGFBR2 mutation analysis


### Introduction

KAPA HiFi DNA polymerase, which was chosen from more than 100 million wild-type DNA polymerase mutants, is a novel PCR enzyme developed and engineered by Kapa Biosystems, exhibiting excellent fidelity and amplification performance.

### Next Generation sequencer (Roche FLX) analysis data

	KAPAHiFi™	P-H from competitor I
Type of enzyme	High fidelity engineered single enzyme system	Mixture of high fidelity enzyme and <i>Taq</i> polymerase
Total sequence number (base)	17,724,794	7,273,424
Total sequence number / Total error number (base)	3,544,959	89,795
Error rate	2.82 E-07	1.11 E-05

\*Data courtesy of Dr. Phillip Buckhaults, University of South Carolina



A huge amount of time and effort is required in order to quantify error rates of high fidelity PCR enzymes. Common methods are “to carry out a great number of sequencing” or “to quantify error rates through blue/white colony selection by cloning E-coli”, both of which are tiring. We performed PCR/sequence analysis of a sequence of about 100bp which is difficult to correctly amplify by PCR, comparing the fidelity of KAPA HiFi DNA Polymerase with competitors’ products in a manner users actually do.

In this Human TGFBR2 mutation analysis, you need to analyze Wild Type with ten consecutive “A” s in a sequence and Mutant with nine consecutive “A” s in a sequence. However, the existence of these consecutive “A” s makes it difficult to perform accurate PCR amplification. In this case, a sequencing using PCR products will not lead to good sequencing results.

In order to verify which PCR reagent brought about the most accurate amplification, we amplified the sequence with a variety of commercially available PCR reagents including KAPA HiFi HotStart DNA Polymerase, and performed sequence analysis.

### Procedure

**Human TGFBR2 mutation analysis was performed in accordance with the following reference literature.**

- 1) Markowitz S, Wang J, Myeroff L, et al. Inactivation of the type II TGF-beta receptor in colon cancer cells with microsatellite instability, *Science* 1995;268:1336-1338.
- 2) Yamashita S et al, Methylation Silencing of Transforming Growth Factor-β Receptor Type II in Rat Prostate Cancers, *Cancer Res* 68:2112-2121, 2008

#### <PCR>

##### ● Primer sequence

Forward AAGCTCCCCTACCATGACT  
 Reverse TGCACATCATCAGAGCTACAG

##### ● Amplification size

118bp

##### ● PCR reagent

- 1 KAPA HiFi HotStart DNA Polymerase
- 2 Ph from competitor F (HotStart)
- 3 Pr from competitor T (HotStart)
- 4 Ko from competitor T (HotStart)
- 5 Wild Type *Taq* DNA Polymerase

##### ● Template DNA

We used the following two types of genomic DNA templates purified from human prostate cancer-derived cell line, whose sequence information is clear.

##### ● Reaction setup

Based on each manufacturer’s recommended conditions

##### ● Thermal cycler

LifePro with gradient function (Bior)

##### ● Cycling program

Based on each manufacturer’s recommended conditions  
 Gradient PCR was performed to determine the optimum annealing temperature.

#### <Sequencing>

Templates for sequencing were prepared by purifying PCR products with a commercially available kit. We performed “Forward” and “Reverse” sequencing of two templates (22Rv1 and LNCaP) respectively.

##### ● Primer sequence

Forward CCCCTACCATGACTTTATTCT  
 Reverse TCATCAGAGCTACAGGAACAC

