

Thermal Stability of Heterotrimeric pMHC Proteins as Determined by Circular Dichroism Spectroscopy

Anna Fuller¹, Aaron Wall¹, Michael D Crowther¹, Angharad Lloyd¹,
Alexei Zhurov², Andrew K. Sewell¹, David K. Cole^{1,*} and Konrad Beck^{2,*}

¹Division of Infection and Immunity and Systems Immunity Research Institute, Cardiff University School of Medicine, Heath Park, Cardiff, UK; ²Cardiff University School of Dentistry, Heath Park, Cardiff, UK

*For correspondence: ColeDK@cf.ac.uk; BeckK@cf.ac.uk

[Abstract] T cell receptor (TCR) recognition of foreign peptide fragments, presented by peptide major histocompatibility complex (pMHC), governs T-cell mediated protection against pathogens and cancer. Many factors govern T-cell sensitivity, including the affinity of the TCR-pMHC interaction and the stability of pMHC on the surface of antigen presenting cells. These factors are particularly relevant for the peptide vaccination field, in which more stable pMHC interactions could enable more effective protection against disease. Here, we discuss a method for the determination of pMHC stability that we have used to investigate HIV immune escape, T-cell sensitivity to cancer antigens and mechanisms leading to autoimmunity.

Keywords: Peptide-MHC stability, Circular dichroism, Thermal stability, T-cells, Peptide vaccines, Recombinant protein, Protein folding

[Background] The ability of CD8⁺ T-cells to respond to foreign invaders or dysregulated self is dependent on stable pMHC class I (pMHCI) presentation at the cell surface. Structurally, MHCI molecules form a peptide binding groove formed of two parallel α helices with a floor of β sheets at the interface between the $\alpha 1$ and $\alpha 2$ domains (Latron *et al.*, 1992). The peptide binding groove has primary peptide binding pockets (B and F) that tightly interact with specific amino acids towards the N- and C-terminals of bound peptides. Although these pockets can accommodate a range of amino acids, they exhibit preferences for certain side chains that have been characterized using structural and biochemical approaches (Parker *et al.*, 1992). This information has been used to generate so called 'heteroclitic' peptides in which natural peptides that have poor MHC-anchors can be modified with amino acids that bind optimally to MHC for vaccination (Cole *et al.*, 2010). Moreover, pMHC stability has been linked to HIV immune escape (Bronke *et al.*, 2013) and the selection of autoreactive T-cell clones (Yin *et al.*, 2011). Thus, understanding the mechanisms that control pMHC stability is important for therapeutic design and understanding complex human diseases. Here, we developed a protocol to accurately determine pMHC stability using circular dichroism spectroscopy. We have used this technique, together with structural, biophysical and cellular experiments, to provide new insight into the molecular factors that determined T-cell antigen recognition in the context of a range of human diseases (Kløverpris *et al.*, 2015; Knight *et al.*, 2015; Motozono *et al.*, 2015; Cole *et al.*, 2016; Jones *et al.*, 2016; Cole *et al.*, 2017).

Materials and Reagents

1. Cellulose nitrate 0.45 µm filter papers (Sartorius, catalog number: 11306-47-N)
2. 1.2 µm glass microfiber filters (GE Healthcare, catalog number: 1822-070)
3. 10 ml plastic syringes, Luer slip BD Plastipak (BD, catalog number: 302188)
4. 25 G needle (BD, catalog number: 300600)
5. 1 ml plastic syringes, Luer slip BD Plastipak (BD, catalog number: 303172)
6. 1.5 ml microcentrifuge tubes
7. Amicon centrifugal concentrating tubes 4 ml MWCO 10 kDa (Merck, catalog number: UFC801096)
8. Phosphate buffered saline (PBS) made up from Dulbecco A tablets (137 mM NaCl, 3 mM KCl, 8 mM Na₂HPO₄, 1.5 mM KH₂PO₄, pH 7.3) (Thermo Fisher Scientific, Thermo Scientific™, catalog number: BR0014G)
9. Ultrapure water (> 18 MΩ cm) for buffer preparations
10. Bolt® Bis-Tris 4-12% precast gels (Thermo Fisher Scientific, Invitrogen™, catalog number: NW04120BOX)
11. BIEye prestained protein markers (Geneflow, catalog number: S6-0024)
12. Quick Coomassie Stain (Generon, catalog number: GEN-QC-STAIN-3L)
13. Ethanol absolute (200 Proof)
14. Nitric acid (HNO₃), 70%

Equipment

1. Reusable bottle top filtration device (Thermo Fisher Scientific, Thermo Scientific™, catalog number: DS0320-5045)
2. Vacuum pump such as KNF Neuberger Vaccum Pump (KNF Neuberger, catalog number: 049268/018121)
3. 500 ml clear Duran bottles (Duran, catalog number: GL 45)
4. Liquid Chromatography system, with a 2 ml injection loop, and a fraction collector; we use the ÄKTA pure 25 L with an F9-R fraction collector (GE Healthcare, model: ÄKTA pure 25 L)
5. Size exclusion chromatography column; we use a Superdex 200 Increase 10/300 GL column, bed volume 24 ml (GE Healthcare, catalog number: 28990944)
6. Benchtop refrigerated micro-centrifuge capable of 14,000 x g (e.g., Eppendorf, model: 5418 R)
7. Far-UV spectrophotometer with a bandwidth < 1.8 nm and quartz cuvettes. We use a single-beam Beckman DU 800 instrument with microcells that allow measurements with volumes of 50 to 100 µl (Beckman Coulter, model: DU® 800)
8. A far-UV circular dichroism (CD) spectrometer with a temperature controlled cell holder. We use an AVIV Model 215 instrument (Aviv Biomedical, model: Aviv Model 215) with a single cell Peltier controlled cell holder, or make use of the Module B end-station spectrophotometer at the B23

Synchrotron Radiation CD (Diamond Light Source, model: B23) Beamline at the Diamond Light Source (Jávorfi *et al.*, 2010; Cole *et al.*, 2016). Alternative instruments are available from Applied Photophysics Ltd (Leatherhead, U.K.), JASCO Inc. (Easton, MD), and Olis Inc. (Bogart, GA)

9. Strain free sealable quartz cuvettes of appropriate path length fitting the CD instrument's cell holder. We use Teflon stoppered Hellma Suprasil cells of various thickness, mostly 0.1-cm
10. A well-ventilated chemical fume hood should be accessible for cleaning cuvettes with HNO₃

Software

1. An analysis software that can import the CD data files and allows curve fitting to a user defined set of equations is required. We use Origin version 7.5 and later (OriginLab Corp., Northampton, MA), but many other programs. *e.g.*, Igor Pro, MATLAB, Micromath Scientist, SigmaPlot, *etc.* will work

Procedure

A. Purification and buffer exchange of pMHC protein

HLA-A2 and human β 2-microglobulin (β 2m) sequences were generated as described and cloned into separate pGMT7 expression plasmids (Cole *et al.*, 2007). Complexes of HLA-A2, β 2m and peptides are refolded and purified as described (Bulek *et al.*, 2012, MacLachlan *et al.*, 2017) with the essential steps summarized below. For CD analysis, proteins should be of high purity and the preparation should lack the presence of denatured or incorrectly folded protein. In order to achieve this, the proteins are purified on the day of, or the day before analysis, and stored on ice to avoid any freeze thaw cycles.

1. Prepare a 500 ml solution of PBS using 1 tablet per 100 ml of ultra-pure water. With the bottle top filtration device and vacuum pump, filter this solution through a cellulose nitrate 0.45 μ m filter paper, with a 1.2 μ m glass microfiber filter as a pre-filter, into a 500 ml Duran bottle. Prepare another 500 ml Duran bottle of ultra-pure water.
2. Connect the gel filtration column to the ÄKTA FPLC system and wash with sterile ultra-pure water followed by equilibrating the system with PBS. Wash the 2 ml injection loop with ultra-pure water and equilibrate with PBS using a 10 ml syringe, taking special care to avoid the occurrence of any air gaps that can be transferred to the column.
3. Using a 25 G needle transfer the protein solution, of no more than 1 ml, into the injection loop using a 1 ml syringe avoiding any air bubbles.
4. Run the column at 0.5 ml/min and monitor elution by UV absorbance at 280 nm (Note 1). Collect 1-ml fractions (the number of fractions depends on the elution profile, but we try and take 2-4 fractions from the center of the peak) and analyse by SDS-PAGE under reducing conditions to check for purity of the pMHC sample. We use 4-12% gradient Bolt® Bis-Tris gels, BLUeye prestained protein ladder, NuPAGE MES running buffer, and Quick Coomassie Stain.

5. Run an SDS-PAGE. Upon SDS-PAGE, pMHC dissociates into its constituent parts. Bands of HLA-A2 and β 2m run at positions corresponding to ~32 and ~12 kDa, respectively. The peptide (~1 kDa) will run out of the gel (Figure 1).
6. Combine only the purest two or three fractions and concentrate using an Amicon centrifugal filtration tube with an MWCO of 10 kDa. This provides protein of adequate purity and concentration for CD analysis.
7. Directly before CD measurements, spin the samples at 14,000 x g for 20 min at 4 °C usually without resulting in any visible pellet.




Figure 1. SDS-PAGE analysis of pMHCI. Approximately 20 μ g of pMHCI (right lane) was loaded onto a 4-12% gradient Bolt® Bis-Tris gel with 3 μ l BLUeye prestained protein ladder (left lane) and run at 200 V for 23 min in NuPAGE MES running buffer. The gel was stained for 5 min in Quick Coomassie Stain and destained with ddH₂O. The pMHCI α chain runs as a separate band at ~32 kDa and β 2m at ~12 kDa.

B. Determination of parameters

1. Record UV spectra of the supernatant (step A7) from 360 to 230 nm in a 1-cm quartz cuvette with the buffer baseline subtracted. Within the 360 to 320 nm range, the absorbance should be essentially zero; any slope is due to light scattering indicating the presence of large aggregates.
2. Calculate the 280 nm absorption coefficient based on the amino acid composition (Pace *et al.*, 1995) using the ProtParam website at <http://web.expasy.org/protparam/>.

Due to the high amount of the aromatic residues Tyr (5.3%) and Trp (3.2%) within the HLA-A2 and β 2-microglobulin sequences combined (see [Supplements 1 and 2](#)), a good estimate for

$A_{280\text{ nm}}^{0.1\%}$ (absorption coefficient) can be expected. As HLA-A2 and β 2-microglobulin have two and one intra-chain disulphide bonds, respectively, the cystine contribution to the absorption is included (ProtParam option: 'assuming all pairs of Cys residues form cystines').

3. Calculate protein concentration from the measured optical density at 280 nm $OD_{280\text{ nm}}^{1\text{ cm}}$ as:

$$c \left(\frac{\text{mg}}{\text{ml}} \right) = \frac{OD_{280\text{ nm}}^{1\text{ cm}}}{A_{280\text{ nm}}^{0.1\%}} \quad (1)$$

4. Use ProtParam to calculate the molecular mass based on the average isotopic masses of amino acids. All three sequences of the HLA-A2/ β 2m/peptide complex are entered into the program input field together. To account for the two missing peptide bonds due to the three instead of single chain complex as seen by the program, two times the molecular mass of water (2 x 18.02) is added to the resulting mass M_r .
5. Calculate the mean residue weight (MRW) required for the normalization of the CD data from this mass by dividing the number of peptide bonds:

$$MRW = \frac{M_r + 2 \times 18.02}{n - 3} \quad (2)$$

where, n is the number of amino acids within the trimer.

C. Data acquisition

1. Based on the concentration of the stock solution determined as described above, dilute an aliquot with buffer to result in *ca.* 600 μ l of sample with $c \sim 0.15$ mg/ml (~ 3 μ M), and record a further absorption spectrum to give the final protein concentration (Note 2).
2. Following the usual start-up procedure of the CD instrument (30 min flushing with N_2 followed by a further 30 min period for stabilizing the xenon arc lamp), record a spectrum of the buffer (step A1) in a clean 0.1-cm stoppered quartz cell from 260 to ~ 195 nm with the cell holder equilibrated to 4 $^{\circ}$ C. This spectrum is termed the baseline. We usually record spectra using a 1 nm bandwidth in 0.2 nm intervals with 3 or 4 sec measuring time per data point resulting in *ca.* 20 to 25 min per spectrum (Note 3).
3. Record a further buffer baseline in the Kinetic Mode at a constant wavelength of 218 nm for 1 min in 2 sec intervals.
4. Clean the cell with water by rinsing, followed by rinsing with absolute ethanol and dried under a stream of nitrogen.
5. Fill the same cell completely with the protein solution (~ 550 μ l), closed tightly with a Teflon stopper avoiding trapping of an air bubble, and place into the CD cell holder in the same orientation as used for recording the baseline.

6. After ~2 min to allow for temperature equilibration, record the protein spectrum using the same instrument settings as used for the baseline. Data are reliable if the photomultiplier dynode voltage is < 500 V (Note 4).
7. For temperature denaturation, set the wavelength to 218 nm (minimum in CD spectrum). Record changes in ellipticities in 0.5 °C intervals from 4 to ~80 °C using the following parameters: data collection 12 sec/point; temperature equilibration 12 sec; temperature dead band 0.3 °C; instrument heating rate 4 °C/min. During the intervals between data collection, the sample is shielded from the light beam by a closed slit. These settings result in an average heating rate of ~36 °C/h (Note 5).
8. At elevated temperatures the HLA-A2/β2m/peptide complexes form large, visible aggregates (precipitation) that result in sharp changes of the melting profile and an increase of the dynode voltage indicating increased scattering. At this point, the measurement can be stopped.
9. Clean protein precipitate containing cuvettes thoroughly rinsing with water followed by filling with 70% HNO₃ in a fume hood, incubation for 2 h to overnight, rinsing with water, ethanol, and drying with a stream of N₂.


Data analysis

1. Import CD data files into the OriginLab software (Note 6).
2. For spectra analysis, subtract the buffer baseline. Normalise θ values recorded in mdeg to molar ellipticities $[\theta]_{MRW}$ in the classical units (deg cm² dmol⁻¹) according:

$$[\theta]_{MRW} = \frac{\theta \times MRW}{c \times d} \quad (3)$$

with the concentration c in mg/ml and cell path length d in mm (Note 7).

3. Subtract the mean of θ recorded at 218 nm in the Kinetic Mode (step C3) from the thermal melt ellipticities, and normalize to $[\theta]_{MRW}^{218nm}$ according to Eq. (3).
4. Analysis of the thermal stability assumes a 2-state mechanism in which three unfolded polypeptide chains U combine to a native complex N :



where, the equilibrium constant K is:

$$K = \frac{c_N}{c_U^3} = \frac{F}{3c_0^2(1-F)^2} \quad (5)$$

with a total chain concentration $c_0 = c_U + 3c_N$ and a degree of conversion to the native state $F =$

$3 c_N/c_0$ (see *e.g.*, Engel *et al.*, 1977; Marky and Breslauer, 1987). The standard Gibbs free energy can be written as:

$$\Delta G^0 = -RT \ln K = \Delta H^0 - T\Delta S^0 \quad (6)$$

in which R is the gas constant ($8.3145 \text{ J K}^{-1} \text{ mol}^{-1}$), ΔH^0 the standard enthalpy, and ΔS^0 the standard entropy. From Eqs. (5) and (6) it follows for the midpoint of transition, where $F = 0.5$ and T equals the melting temperature T_m , at which half of the proteins are in the folded and unfolded state, respectively, that

$$K = 1/0.75 c_0^2 \text{ and } T_m = \frac{\Delta H^0}{\Delta S^0 + R \ln(0.75 c_0^2)} \quad (7)$$

This means that in contrast to a single chain polypeptide the melting temperature is concentration dependent. T_m and ΔH^0 can be obtained by fitting the entire transition curve to the van't Hoff equation:

$$K = \exp \left[\frac{\Delta H^0}{RT} \left(\frac{T}{T_m} - 1 \right) - \ln 0.75 c_0^2 \right] \quad (8)$$

obtained by solving Eq. (7) for ΔS^0 and substituting for ΔS^0 in Eq. (6).

The molar ellipticities for the native and unfolded state $[\Theta]_n$ and $[\Theta]_u$ show an additional temperature dependence which in a first approximation can be assumed as linear:

$$[\Theta]_n = [\Theta]_n^0 + b_n T \text{ and } [\Theta]_u = [\Theta]_u^0 + b_u T \quad (9)$$

where, $[\Theta]_n^0$ and $[\Theta]_u^0$ correspond to the $[\Theta]$ intercept at $T = 0$, and b_n and b_u describe the slopes of the native and unfolded states, respectively. The entire transition curve is thus described by:

$$[\Theta] = F ([\Theta]_n^0 + b_n T - [\Theta]_u^0 - b_u T) + [\Theta]_u^0 + b_u T \quad (10)$$

which is used for fitting with the Levenberg-Marquardt algorithm implemented in OriginLab. The essential parts of the fitting script are included as [Supplement 4](#).

- As the HLA-A2/β2m/peptide complexes tend to precipitate within the transition region (step C8), we assume a common molar ellipticity $[\Theta]_u^0 = -4,500 \text{ deg cm}^2 \text{ dmol}^{-1}$ and slope $b_u = 0$. For many temperature denatured proteins and model peptides $[\Theta]_u^0$ values in the $-4,000$ to $-5,000 \text{ deg cm}^2 \text{ dmol}^{-1}$ range are observed (Venyaninov *et al.*, 1993) (Notes 8 and 9).
- Representative data

The general procedure to analyse the stability of a trimeric pMHC complex is illustrated using HLA-A2/β2m/ILAKFLHWL as an example (Figure 2). This complex has been analysed in detail

(Cole *et al.*, 2017) and the structure has been solved in complex with a TCR (Protein Data Bank entry 5MEN).




Figure 2. Thermal stability of HLA A2/β2m/ILKA. A. Protein concentration is determined by UV absorbance spectroscopy. Wavelength scans from 360 to 230 nm are recorded in a 1-cm quartz cuvette using the same sample and buffer as used in subsequent CD measurements. Concentration is calculated assuming an absorbance coefficient $A_{280\text{ nm}}^{0.1\%} = 2.257$ as calculated based on the amino acid composition. B. A far-UV CD spectrum was recorded in a 0.1-cm quartz cuvette at 4 °C and buffer baselines measured in the same cell using the same instrument parameters were subtracted. The signals were normalized to mean-residue-weight (MRW) ellipticities $[\theta]_{MRW}$ acc. Eq. (3). C. Thermal stability was determined by following the ellipticity Θ at 218 nm upon increasing the temperature with measuring 12 sec per data point in $\Delta T = 0.5$ °C intervals, an equilibration time of 12 sec between data points, a temperature dead band of ± 0.3 °C, and using slit closure between measurements resulting in a mean heating rate of 36 °C/h (red, left axis). The sample forms large, visible aggregates indicating precipitation in the 58 to 65 °C range (open circles). Data points below this temperature (filled circles) were fitted acc. Eq. (10) resulting in the dashed curve. Aggregate formation is accompanied by increasing light scattering resulting in a steep increase of the dynode voltage (blue, right axis) followed by a decrease when large particles sink to the bottom of the cuvette.

- a. Absorption coefficient and mean-residue-weight are calculated by putting all three sequences ([Supplements 1 to 3](#)) into the input frame of ProtParam (step B2) resulting in:

Number of amino acids: 386

Molecular weight: 44820.03

.....

Extinction coefficients are in units of $M^{-1} cm^{-1}$, at 280 nm measured in water.

Ext. coefficient 101675

Abs 0.1% (=1 g/L) 2.269, assuming all pairs of Cys residues form cystines

- b. The measured absorption of the protein solution (Figure 2A) shows $OD_{280 nm}^{1 cm} = 0.283$ resulting with Eq. (1) in a concentration of $0.283/2.269 = 0.125$ mg/ml.
- c. MRW is calculated according to Eq. (2) as $MRW = (44,820.03 + 2 \times 18.02)/(386 - 3) = 117.12$.
- d. A CD buffer baseline and spectra of both solutions are recorded according to steps C2 to C6 using a 0.1-cm path length cuvette. Baselines were subtracted, and the spectrum is normalised according to Eq. (3) (Figure 2B). At 195 nm, the dynode voltage was ~500 V. This reflects the high absorbance of the chloride ion included in PBS. The spectrum shows a minimum of $[\theta]_{MRW} = -8,950$ deg cm^2 $dmol^{-1}$ at 218 nm.
- e. The melting curve is recorded as described in step C7 from 4 to 65 °C, the mean value of the baseline (step C3) is subtracted, and molar residue ellipticities calculated according to Eq. (3) (Figure 2C). Using the OriginLab script given as [Supplement 4](#) fit the data according to Eq. (10). Values in the 59 to 65 °C range that show a pronounced deviation from a homogenous quasi-sigmoidal curve indicating protein aggregation (open dots in Figure 2C) were omitted from the fits. Fitting results in $T_m = 54.5 \pm 0.1$ °C and $\Delta H^0 = -500 \pm 20$ kJ/mol with a coefficient of determination $r^2 = 0.988$.

Notes

1. Peak fractionation of the pMHC will occur at an elution volume of ~12 ml; it is important to keep only the purest fractions and discard any that contain protein from more than one peak.
2. To determine thermodynamic parameters, for HLA-A2/β2m/peptide complexes we aim for a concentration of ~0.15 mg/ml as for most peptides studied so far this can be relatively easily achieved and provides for a good signal-to-noise ratio to follow CD transition curves at 218 nm. Due to the concentration dependence of T_m for multi-chain proteins [Eq. (7)], it is important to settle for a narrow range to allow for comparison of data.
3. A cuvette can be regarded as optically clean if the baseline using PBS or water runs about parallel to a spectrum recorded without a cell in the light beam ('air baseline'); minima observed

in the 200 to 225 nm range usually indicate that protein from previous measurements is absorbed to the quartz surface.

4. On instruments of other manufacturers, the dynode voltage is called HT (high tension), and the reliable voltage range can differ depending on the type of photomultiplier and electronics.
5. The heating parameters have been tested to result in the set equilibrium temperatures by measuring with a Pt100 Resistance Temperature Detector inside a 0.1-cm cuvette.
6. The AVIV instrument software stores data in ASCII format. Some instruments from other manufacturers use propriety formats that require exporting into ASCII for analysis via other programs.
7. $[\theta]_{MRW}$ can be converted into molar absorption units $\Delta\epsilon = \epsilon_L - \epsilon_R = [\theta]_{MRW}/3,298$, where $\Delta\epsilon$ is the absorption difference for left- and right-handed circularly polarized light.
8. From the various MHC complexes studied in our lab, HLA-A*0101/ β 2m/peptide VTEHDTLLY exhibit high stability both with respect to T_m (68.9 °C) and enthalpy (-719 kJ mol⁻¹); precipitation occurred only at $T > 72$ °C, $[\Theta]_U^0 > -5,200$ deg cm² dmol⁻¹ (Jones *et al.*, 2016). Using $[\Theta]_U^0 = -4,500$ deg cm² dmol⁻¹ resulted in a $r^2 > 0.996$.
9. For strict thermodynamic analysis reversibility of the transition [Eq. (4)] is a prime requirement. It has been suggested, however, that even if reversibility is limited analysis is still of value as temperature induced cooperative unfolding is usually a much faster process than the aggregation of unfolded protein (Privalov, 2009). To allow for comparison of data, we use the same heating parameters for all samples (see step C7).

Acknowledgments

DKC is a Wellcome Trust Research Career Development Fellow (WT095767). AKS is a Wellcome Trust Senior Investigator. This and original work using the described procedures was supported by the UK Biotechnology and Biological Sciences Research Council (BBSRC) (Grant BB/H001085/1). Purchase of the CD instrument was partially funded by BBSRC grant 75/REI18433.

References

1. Bronke, C., Almeida, C. A., McKinnon, E., Roberts, S. G., Keane, N. M., Chopra, A., Carlson, J. M., Heckerman, D., Mallal, S. and John, M. (2013). [HIV escape mutations occur preferentially at HLA-binding sites of CD8 T-cell epitopes.](#) *AIDS* 27(6): 899-905.
2. Bulek, A. M., Cole, D. K., Skowera, A., Dolton, G., Gras, S., Madura, F., Fuller, A., Miles, J. J., Gostick, E., Price, D. A., Drijfhout, J. W., Knight, R. R., Huang, G. C., Lissin, N., Molloy, P. E., Wooldridge, L., Jakobsen, B. K., Rossjohn, J., Peakman, M., Rizkallah, P. J. and Sewell, A. K. (2012). [Structural basis for the killing of human beta cells by CD8⁺ T cells in type 1 diabetes.](#) *Nat Immunol* 13(3): 283-289.

3. Cole, D. K., Bulek, A. M., Dolton, G., Schauenberg, A. J., Szomolay, B., Rittase, W., Trimby, A., Jothikumar, P., Fuller, A., Skowera, A., Rossjohn, J., Zhu, C., Miles, J. J., Peakman, M., Wooldridge, L., Rizkallah, P. J. and Sewell, A. K. (2016). [Hotspot autoimmune T cell receptor binding underlies pathogen and insulin peptide cross-reactivity](#). *J Clin Invest* 126(6): 2191-2204.
4. Cole, D. K., Edwards, E. S., Wynn, K. K., Clement, M., Miles, J. J., Ladell, K., Ekeruche, J., Gostick, E., Adams, K. J., Skowera, A., Peakman, M., Wooldridge, L., Price, D. A. and Sewell, A. K. (2010). [Modification of MHC anchor residues generates heteroclitic peptides that alter TCR binding and T cell recognition](#). *J Immunol* 185(4): 2600-2610.
5. Cole, D. K., Pumphrey, N. J., Boulter, J. M., Sami, M., Bell, J. I., Gostick, E., Price, D. A., Gao, G. F., Sewell, A. K. and Jakobsen, B. K. (2007). [Human TCR-binding affinity is governed by MHC class restriction](#). *J Immunol* 178(9): 5727-5734.
6. Cole, D. K., van den Berg, H. A., Lloyd, A., Crowther, M. D., Beck, K., Ekeruche-Makinde, J., Miles, J. J., Bulek, A. M., Dolton, G., Schauenburg, A. J., Wall, A., Fuller, A., Clement, M., Laugel, B., Rizkallah, P. J., Wooldridge, L. and Sewell, A. K. (2017). [Structural mechanism underpinning cross-reactivity of a CD8⁺ T-cell clone that recognizes a peptide derived from human telomerase reverse transcriptase](#). *J Biol Chem* 292(3): 802-813.
7. Engel, J., Chen, H. T., Prockop, D. J. and Klump, H. (1977). [The triple helix in equilibrium with coil conversion of collagen-like polytripeptides in aqueous and nonaqueous solvents. Comparison of the thermodynamic parameters and the binding of water to \(L-Pro-L-Pro-Gly\)_n and \(L-Pro-L-Hyp-Gly\)_n](#). *Biopolymers* 16(3): 601-622.
8. Jávorfí, T., Hussain, R., Myatt, D. and Siligardi, G. (2010). [Measuring circular dichroism in a capillary cell using the b23 synchrotron radiation CD beamline at diamond light source](#). *Chirality* 22 Suppl 1: E149-153.
9. Jones, K., Wockner, L., Brennan, R. M., Keane, C., Chattopadhyay, P. K., Roederer, M., Price, D. A., Cole, D. K., Hassan, B., Beck, K., Gottlieb, D., Ritchie, D. S., Seymour, J. F., Vari, F., Crooks, P., Burrows, S. R. and Gandhi, M. K. (2016). [The impact of HLA class I and EBV latency-II antigen-specific CD8⁺ T cells on the pathogenesis of EBV⁺ Hodgkin lymphoma](#). *Clin Exp Immunol* 183(2): 206-220.
10. Kløverpris, H. N., Cole, D. K., Fuller, A., Carlson, J., Beck, K., Schauenburg, A. J., Rizkallah, P. J., Buus, S., Sewell, A. K. and Goulder, P. (2015). [A molecular switch in immunodominant HIV-1-specific CD8 T-cell epitopes shapes differential HLA-restricted escape](#). *Retrovirology* 12: 20.
11. Knight, R. R., Dolton, G., Kronenberg-Versteeg, D., Eichmann, M., Zhao, M., Huang, G. C., Beck, K., Cole, D. K., Sewell, A. K., Skowera, A. and Peakman, M. (2015). [A distinct immunogenic region of glutamic acid decarboxylase 65 is naturally processed and presented by human islet cells to cytotoxic CD8 T cells](#). *Clin Exp Immunol* 179(1): 100-107.
12. Latron, F., Pazmany, L., Morrison, J., Moots, R., Saper, M. A., McMichael, A. and Strominger, J. L. (1992). [A critical role for conserved residues in the cleft of HLA-A2 in presentation of a nonapeptide to T cells](#). *Science* 257(5072): 964-967.

13. MacLachlan, B.J., Greenshields-Watson, A., Mason, G.H., Schauenburg, A.J., Bianchi, V., Rizkallah, P.J., Sewell, A.K., Fuller A. and Cole, D.K. (2017). [Using X-ray crystallography, biophysics, and functional assays to determine the mechanisms governing T-cell receptor recognition of cancer antigens.](#) *J Vis Exp.* 6(120).
14. Marky, L. A. and Breslauer, K. J. (1987). [Calculating thermodynamic data for transitions of any molecularity from equilibrium melting curves.](#) *Biopolymers* 26(9): 1601-1620.
15. Motozono, C., Pearson, J. A., De Leenheer, E., Rizkallah, P. J., Beck, K., Trimby, A., Sewell, A. K., Wong, F. S. and Cole, D. K. (2015). [Distortion of the major histocompatibility complex class I binding groove to accommodate an insulin-derived 10-Mer peptide.](#) *J Biol Chem* 290(31): 18924-18933.
16. Pace, C. N., Vajdos, F., Fee, L., Grimsley, G., and Gray, T. (1995). [How to measure and predict the molar absorption coefficient of a protein.](#) *Protein Sci* 4(11): 2411-2423.
17. Parker, K. C., Bednarek, M. A., Hull, L. K., Utz, U., Cunningham, B., Zweerink, H. J., Biddison, W. E. and Coligan, J. E. (1992). [Sequence motifs important for peptide binding to the human MHC class I molecule, HLA-A2.](#) *J Immunol* 149(11): 3580-3587.
18. Privalov, P. L. (2009). [Microcalorimetry of proteins and their complexes.](#) *Methods Mol Biol* 490: 1-39.
19. Venyaminov S. Yu., Baikalov I. A., Shen, Z. M., Wu, C. S., and Yang, J. T. (1993). [Circular dichroic analysis of denatured proteins: inclusion of denatured proteins in the reference set.](#) *Anal Biochem.* 214(1):17-24.
20. Yin, Y., Li, Y., Kerzic, M. C., Martin, R. and Mariuzza, R. A. (2011). [Structure of a TCR with high affinity for self-antigen reveals basis for escape from negative selection.](#) *EMBO J* 30(6): 1137-1148.